

BMC ENTRY  
LEVEL MODEL

# FIXED PRESSURE CPAP MACHINE

CPAP Machines provide a constant (fixed) pressure for the clinical treatment of Obstructive Sleep Apnea. They need to be set up to your specific pressure requirements, based on your Doctor's prescription or existing equipment machine configuration. Typically CPAP machines start off at a low initial pressure which then slowly increases to the specified treatment pressure over a predetermined (Ramp) time, which is roughly the time it takes to fall asleep. Once the machine reaches the treatment pressure it remains there for the remainder of the sleep period.



## Luna iQ – CPAP

- 🔗 Latest BMC Model
- 🔗 Integrated Humidifier
- 🔗 Quietest ever machine from BMC
- 🔗 Fixed Pressure
- 🔗 Therapy Software Included on SD Card
- 🔗 Remote Sleep Reports Included
- 🔗 5 Year Machine Warranty
- 🔗 Size : 26x17x10cm (LxWxH)
- 🔗 Weight : 1.8kg

## ADD A BMC, PHILIPS OR F&P MASK



- 🔗 BMC P2 Nasal Pillow Mask
- 🔗 BMC N5 Micro Nasal Mask
- 🔗 BMC F1A Full Face Mask

- 🔗 F&P Brevida Nasal Pillow Mask
- 🔗 Philips DreamWear UtN and Gel Pillow Nasal Masks
- 🔗 F&P Simplus Full Face Mask



[cpapsales.com.au](http://cpapsales.com.au) | [support@cpapsales.com.au](mailto:support@cpapsales.com.au) | 1300 04 2727

*This material is provided for informational and promotional purposes and should not be considered as a substitute for professional medical advice, diagnosis or treatment.*